

# cpaphero



cpap cushion cleanser



- **Included:** 3 month supply of cpaptabs - 96 count
- **Benefits:** Gently cleanses and reduces buildup of dirt, debris, grime, and more on your cpap cushions.
- **Compatibility:** Works with all cushions for cpap masks that are directed to be submerged in water to clean.
- Use with **nasal pillows, nasal or full face** cpap cushions
- Simply fill a glass with warm water, place the tab in the water, place the cushion in the water. Wait 3-5 minutes and watch the water change colors. Remove and rinse and they are ready to use