

RhythmStar Battery Usage & Charging

Each RhythmStar Mobile Cardiac Monitor is supplied with two batteries and a charger. To ensure you always have power available we recommend charging one battery whenever the other is in use. A fresh battery will generally last 24-hours; and a drained battery takes about 4-hours to charge.

Check your monitor's battery level indicator regularly; and change batteries as needed at least once a day. Certain conditions, like poor cellular network coverage or high data usage, may increase battery drain and require more frequent battery changes.

To charge a battery, place it in the charger as shown, and plug the charger into a wall outlet. When a battery is in place, the charger's indicator light will glow **solid red** to indicate charging or **solid blue** to indicate a full charge. (The charger will NOT overcharge a battery.)



**Red
Charging**



**Blue
Fully Charged**